

URIT

A better choice for portable & reliable
urinalysis with multifunction

URIT-31 Urine Analyzer



Global Diagnostics Supplier

URIT-31

Urine Analyzer

Specifications

Parameters: pH, Nitrite, Specific Gravity, Blood, Glucose, Bilirubin, Urobilinogen, Ketone, Leucocytes, Protein, Ascorbic Acid, Creatinine, Microalbumin, Calcium

Measurement Principle: dual wavelength reflectance method

Measuring Wavelength: 4 LED wavelengths (470, 550, 620, 720nm)

Sample preparation: manual dipping

Reaction time: 60 seconds per test

Data memory: up to 1000 samples

Display: LCD colored touch-screen

Strips: URIT 10G/11G (14G optional) urine reagent strips

Interface: Micro USB, Bluetooth

Work condition: environmental temperature 5–40°C, RH:20%~80%

Power supply: power adapter 9V DC, 1.5A

Dimensions: 120(L) × 83(W) × 31(H)mm (analyzer)
110(L) × 80(W) × 38(H)mm (printer)

Weight: 180g (analyzer), 180g (printer)



Features

- Point-of-Care urinalysis with palm-sized feather weight units
- Bluetooth technology for wireless printing and data transfer
- LCD colored touch-screen
- Single button easy operation
- Recharge, transfer data and upgrade, all with Micro USB cable
- Compatible with URIT 10G, 11G and 14G urine reagent strips



URIT URINE CONTROL



URIT 10G,11G,14G urine reagent strip

14G optional with Microalbumin, Creatinine & Calcium

URIT MEDICAL ELECTRONIC (GROUP) CO., LTD.

Head Office

Add.: No.4 East Alley, Jiu Hua Road, Guilin,
Guangxi 541001, P.R. China

Tel: +86-773-2288586 Fax: +86-773-2288560

E-mail: export@uritest.com service@uritest.com

Shenzhen Office

Add.: Room E.J, 24th Floor, New Baohui Building, Nanshan Road,
Nanshan District, Shenzhen, P.R. China

Tel: +86-755-26050853 Fax: +86-755-86036150

Http://www.urit.com



ISO9001
ISO13485